

Programmable and customizable assistance

Electronics by Triride

USB LCD Display with 5 levels of assistance

Double brake (Disk brake e V-brake)

BLUETOOTH SYSTEM: IBS - ACC - Reverse Advance Battery Regen

Start & Walk Assist

¹Weight refers to device excluding battery and arms for wheelchair connection.

² Autonomy refers to 36 V battery on a flat route at medium speed of 6 km/h (level 1) with user weight 70 kg. The autonomy depends on user weight, road conditions, battery model, speed/driving mode and on the general state of the wheelchair and tires.

Triride srl reserves the right to modify the models with updates and innovations and therefore the images may differ from the final product sold.